




TRIMETHOXSILANES FOR SILICA GEL BONDING

NEW!



UCT SPECIALTIES, LLC • 2731 Bartram Road, Bristol, PA 19007
Ph: 215-781-9255 • Fax: 215-785-1226 • www.unitedchem.com

A0698

(Aminoethylaminomethyl)phenethyltrimethoxysilane

QTY: 3.00 KG

LOT #: X

CAS#: 74113-77-2




HEALTH	3
FLAMMABILITY	1
REACTIVITY	1

FOR CHEMICAL PURPOSES ONLY • NOT FOR DRUG USE

PETRARCH® Chemicals, manufactured by UCT Specialties announces the availability of several unique trimethoxysilanes. These materials can be bonded to the surface of silica gel particles in a similar manner to other methoxy silanes that are used to create traditional C18, cyano, or other phases seen in HPLC. Once bonded, these phases offer analysts unique chemistries for the retention and separation of a range of challenging compounds.

INNOVATION THROUGH CHEMISTRY


N-(3-(Trimethoxysilyl)propyl)dimethylammoniopropylsulfonate



Zwitterion – may be useful and show interesting selectivity for mixtures of acids, bases, and neutral compounds.


Part Number	CAS Number	Formula	MW	BP/mm/Hg	SG	nd
T2912	N/A	C ₁₁ H ₂₇ NO ₆ SSi	329.2	N/A	N/A	N/A

(Aminoethylaminomethyl)phenethyltrimethoxysilane




Part Number	CAS Number	Formula	MW	BP/mm/Hg	SG	nd
A0698	74113-77-2	C ₁₄ H ₂₆ N ₂ O ₃ Si	298.46	126-130/0.2	1.02	1.5083

3-Piperazinopropyltrimethoxysilane



Part Number	CAS Number	Formula	MW	BP/mm/Hg	SG	nd
P0686	40762-28-5	C ₁₀ H ₂₄ N ₂ Si	248.39	108°C/1	0.985	N/A

3-Piperidinopropyltrimethoxysilane



Part Number	CAS Number	Formula	MW	BP/mm/Hg	SG	nd
P0683	76750-52-2	C ₁₁ H ₂₅ O ₃ NSi	247.41	N/A	N/A	N/A

Piperazine and piperidine based compounds are commonly found throughout the pharmaceutical and nutraceutical industries.

PRICES AND TERMS

Our prices are subject to change without notice. The price in effect when we receive your order will apply. All prices are in US Dollars and are F.O.B. Bristol, PA 19007. Terms of payment are net 30 days.

MINIMUM ORDERS

We welcome all orders, therefore, we do not have a minimum order requirement. When ordering, please include your purchase order number, complete "Ship To" and "Bill To" address, catalog number, quantity, and description of product(s). Also include your name and a phone number where you can be reached should we have any questions concerning your order.

SHIPMENTS

Normal process is shipment within one business day after receipt of an order if material is in stock. Unless special shipping arrangements are made, all orders will be sent via Federal Express Ground. Shipping and handling charges (freight, insurance, hazardous materials fee) will be added to the invoice unless requested by the customer.

SPECIAL PRICING

We offer special pricing for volume purchases and standing orders. Please call a sales representative for more information on special pricing qualifications.

RETURN POLICY

Our sales staff will process your return. Before returning merchandise, please call to obtain a return authorization number from a sales representative. We will need to know the reason for the return, date of purchase, purchase order number, invoice number, catalog number and lot number of product returned in order to issue a return authorization number.

WARRANTY

All products manufactured by UCT are guaranteed against defects in materials and workmanship for a period of 90 days after shipment. UCT will replace any items that prove to be defective during this time period. The exclusive remedy requires the end user to first advise UCT of the defective product by phone or in writing and must include order number, the lot number and the shipping date.

To initiate this action, photographs of the product, including packaging and labeling of the containers, must be submitted to the UCT Representative for approval. With approval a return authorization can be initiated, and must be received within 30 days. Once the materials arrive at UCT a further inspection of the materials must be completed and accepted by our Quality Manager prior to further action of credits or replacement. UCT's total liability is limited to the replacement cost of UCT products.

This warranty does not apply to damage resulting from misuse.

Contact Us

Phone: 215.781.9255
800.541.0559
Fax: 215.785.1226

UCT, Inc.
2731 Bartram Rd.
Bristol, PA 19007

Email: uct_specialties@unitedchem.com
Web: www.unitedchem.com

